

See discussions, stats, and author profiles for this publication at: <https://www.researchgate.net/publication/228907413>

Miniaturized FTIR-spectrometer Based on Optical MEMS Translatory Actuator

Article in *Proceedings of SPIE - The International Society for Optical Engineering* · January 2007

DOI: 10.1117/12.697898

CITATIONS

37

READS

1,213

5 authors, including:



Thilo Sandner

Fraunhofer Institute for Photonic Microsystems IPMS

182 PUBLICATIONS 1,039 CITATIONS

[SEE PROFILE](#)



Christian Drabe

Fraunhofer Institute for Photonic Microsystems IPMS

33 PUBLICATIONS 292 CITATIONS

[SEE PROFILE](#)



W. Scherf

Carinthian Tech Research AG

26 PUBLICATIONS 518 CITATIONS

[SEE PROFILE](#)




Some of the authors of this publication are also working on these related projects:



Micro Scanner [View project](#)



TACO: Three-dimensional Camera with Object Detection and Foveation [View project](#)



Miniaturized FTIR-Spectrometer based on Optical MEMS Translatory Actuator

Thilo Sandner, Christian Drabe, Harald Schenk

Andreas Kenda*, Werner Scherf*

Fraunhofer Institute for Photonic Microsystems, Dresden, Germany

**Carinthian Tech Research AG, Villach, Austria*

I Introduction

- Advantages of FT spectrometers
- Motivation for MOEMS based FT spectrometers

II MEMS optical path length modulator

- Translatory mirror design
- Experimental results

III System integration


- FTIR System description
- Interferogram sampling concept
- Results / sample spectra

IV Summary

MEMS Scanning Mirror

Properties


- Large optical scan range up to 112°
- High frequencies up to 32 kHz
- Mirror diameter: 0.5 - 3 mm
- 1D and 2D deflection
- High vibration and shock stability
- Excellent long-run behaviour
- Qualified fabrication process
- resonant operation



■ Introduction to FT spectrometry

4

- The main component of a FT spectrometer is a two beam interferometer e.g. a Michelson interferometer.
- The output signal of monochromatic radiation is a cosine wave or of a polychromatic source it is the sum of all cosine waves which is called interferogram.
- A numerical cosine Fourier transformation from time domain (interferogram) to frequency domain yields the spectrum.
- The spectral resolution $\Delta\nu$ depends on the maximum possible OPD of the interferometer.



Schematics of a Michelson interferometer

$$\Delta\nu \approx \frac{1}{2OCD}$$

■ FT versus diffractive spectrometry

5

Advantages

- *Fellgett advantage* (multiplex advantage): All frequencies are measured simultaneously.
- *Jaquinot advantage*: No optical stops -> higher energy throughput.
- *Connes advantage*: Laser referenced wavelength classification.
- Wide spectral range accessible: Limited by dispersion and detector response only.
- Constant resolution over the entire spectral range.

Drawbacks

- Complex, high precision mirror drives.
- Shock sensitive instrumentation.



Development objective / motivation

6

Advantages due to MOEMS technology

Replacement of the macroscopic mirror and its drive with an oscillating micro-mirror.


- Increased reliability and ruggedness.
- System miniaturization.
- Cost reduction.
- Ultra-rapid scan capability. Acquisition time of 0.2 ms for a single scan.

Example of MEMS based FTS:


a) Translatory mirror, b) Lamellar grating booth with in-plane-comb drive

Quelle:

O. Manzardo, Ph.D. Thesis, Neuchatel, 2002



a)



b)

I Introduction

- Advantages of FT spectrometers
- Motivation for MOEMS based FT spectrometers

II MEMS optical path length modulator

- Translatory mirror design
- Experimental results


III System integration

- FTIR System description
- Interferogram sampling concept
- Results / sample spectra



IV Summary

Translatory MOEMS for FT-IR


8



- Fabricated with bulk micromachining technology using 100μm SOI.
- Folded springs to reduce mechanical stress.
- Use for optical path length modulation.



MEMS-Process ame1




Comb-Drive for resonant operation

10

How does it work?


- Electro-statically driven by rectangular-shaped voltage applied to in-plane inter-digital comb fingers.
- Sense cross-over (optical)
- Switch off drive voltage for $\frac{1}{4}$ period

→ large deflection at low voltage




Translatory MOEMS for FT-IR

11




- Mech. amplitude +/- 100 μm ⇒ resolution > 25 cm⁻¹ (single sided acquisition)
- Oscillation frequency ~5 kHz ⇒ interferogram sampling bandwidth 1...10 MHz
- Mirror area ~1.6 mm²
- Mirror deformation 200 nm (RMS)
- Driving voltage 36 V

- At low pressures (10 Pa - 100 Pa) the oscillation starts with an increasing frequency sweep.
- The amplitude grows with increasing frequency.
- The oscillation breaks down as the mech. oscillation gets in phase with the driving signal.
- Due to the strong damping no oscillation was detected at higher pressure levels.




Amplitude versus frequency at 10 Pa




Maximum amplitude versus voltage for varying pressure

■ Novel translatory MOEMS for FT-IR

13



MOEMS mirror device (new design)



Novel design of mirror suspension (Pantograph)

■ Objectives for novel MOEMS device


Increased deflection (OPD) up to 1 mm $\Rightarrow \Delta\nu=10\text{cm}^{-1}$

Increased aperture up to 7 mm² (will raise SNR)


Reduced oscillation frequency (500 Hz) to reduce electronical bandwidth

■ Novel translatory MOEMS for FT-IR

14



SEM of Pantograph suspension




Detail

Parameters of novel translatory MOEMS device


- Large mirror aperture of $\varnothing = 3 \text{ mm}$
- Resonant oscillation @ 500 Hz / 1000 Hz
- Target amplitude up to $\pm 250 \text{ }\mu\text{m}$ @ pressures $< 500 \text{ Pa}$
- Fabrication in standard IPMS scanner process using $30 \text{ }\mu\text{m}$ SOI

Preliminary results:

Amplitude versus Frequency (1 kHz)



Amplitude versus Voltage (1 kHz)




Measurement Setup


- Michelson interferometer
- Ambient pressure: 10 Pa - 500 Pa
- Driving voltage: 10 V - 42 V

Experimental Results

Amplitude versus Frequency @ 1kHz, 20Pa



Amplitude versus Druck @ 500Hz



- Maximale Amplitude von $\pm 240\mu\text{m}$ bzw. OPD von $960\mu\text{m}$ @ 20Pa & 100V
- Deutliche Reduzierung der Squeeze-Filmdämpfung durch erhöhte Rückseitenkavität
 $\rightarrow \pm 95\mu\text{m}$ Amplitude @ 10kPa (40V) bzw. $\pm 80\mu\text{m}$ @ 100kPa
- Optimierung der Modentrennung erforderlich

I Introduction

- Advantages of FT spectrometers
- Motivation for MOEMS based FT spectrometers

II MEMS optical path length modulator

- Translatory mirror design
- Experimental results


III System integration

- FTIR System description
- Interferogram sampling concept
- Results / sample spectra


IV Summary

MOEMS FTIR system schematics

18



Optical layout and block diagram of the signal path of the FTIR system.





- Fibre coupled optics (Chalcogenide fibre).
- Operation of the MOEMS device within a vacuum chamber.
- Detector: Thermoelectrically cooled photovoltaic MCT.
- CaF_2 Beamsplitter and windows.
- Reference Laser interferometer accesses the backside of the MOEMS device.
- Laser source: Temp. Stabilized VCSEL $\lambda=760\text{nm}$.
- DSP based data acquisition and processing.
- Dispersive optical elements and detector response define wavelength range to $\sim 2\text{...}6\text{ }\mu\text{m}$.



- # Fraunhofer

MOEMS FTIR prototype

21



Fraunhofer




thilo.sandner@ipms.fraunhofer.de

Spectrometer Workshop
Jena, 2008-03-13




Interferogram sampling

22



Reference signal clocked sampling eliminates distortions due to the mirror's nonlinear motion.


Interferogram of polychromatic infrared radiation



Spectrum (FFT Power plot)

24

Staggered plot of infrared spectra with and without averaging




Averaging of single scans
raises SNR a $(n)^{0.5}$

- SNR 240, 24dB
- SNR 8 , 9dB

Transmission plot

1.5mm Polystyrene



Absorbance plots

