


Quantitative Laser Diagnostics for Combustion Chemistry and Propulsion

Lecture 9: Tunable Diode Laser Absorption (TDLAS)

1. History and vision – for aeropropulsion
2. Absorption fundamentals
3. Absorption sensor strategies
4. Wavelength access – lasers and detectors
5. Example applications – aerospace
6. Future trends for aerospace

Direct-connect scramjet combustor
at UVa flow facility



1. The History of TDL Absorption for Aeropropulsion: 35 Years: From the Laboratory to Flight

- 1977 – TDL absorption in shock tube flows and flames
- 1989 – Mass flux sensor using O_2 absorption
- 1993 – Multiplexed measurements of H_2O , T and momentum flux
- 1998 – Combustion control (lab flames, incinerator)
- 2001 – Multi-species in flames: CO , CO_2 , NH_3 , H_2O
- 1996-present – Applications to flow facility characterization: arcjets, hypersonic flow tunnels, gas turbine engine sector rigs...
- 1998-present – Applications for engine tests: scramjet combustors, commercial aircraft engines, ic-engines, pulse detonation engines, gas turbines, augmentors...
- 2012 – TDL absorption by AFRL for T and combustion efficiency in scramjet flight test

T/species in Shocktube




Hanson, *Appl Opt* (1977)

SCRAMJET @ WPAFB



Rieker, *Proc Comb Inst* (2009)


HiFire-2 Scramjet Flight Test



AFRL WPAFB: scheduled 2012

1. Vision for TDLAS Sensors for Aeropropulsion

- Diode laser absorption sensors offer prospects for time-resolved, multi-parameter, multi-location sensing for performance testing, model validation, feedback control





- Sensors developed for T, V, H_2O , CO_2 , O_2 , & other species
- Prototypes tested and validated at Stanford
- Several successful demonstrations in ground test facilities
- Opportunities emerging for use in flight

Now for some absorption fundamentals



2. Absorption Fundamentals: The Basics

Absorption of monochromatic light





- Scanned-wavelength *line-of-sight* direct absorption
 - Beer-Lambert relation $\tau_\nu \equiv \frac{I_t}{I_o} = \exp(\underbrace{-k_\nu \cdot L}_{\text{absorbance}}) = \exp(-n_i \cdot \sigma_\nu \cdot L)$
 - Spectral absorption coefficient $k_\nu = S(T) \cdot \Phi(T, P, \chi_i) \cdot \chi_i \cdot P$

2. Absorption Fundamentals: The Basics





- Shifts & shape of Φ contain information (T, V, P, χ_i)

2. Absorption Fundamentals: The Basics



- T from ratio of absorption at two wavelengths

2. Absorption Fundamentals: Summary




- Wavelength multiplexing is also effective
 - To monitor multiple parameters or species
 - To assess non-uniformity along line-of-sight

Two primary strategies for absorption measurements

3. Absorption Sensor Strategies:


Direct Absorption (DA) & Wavelength Modulation Spectroscopy (WMS)



- Direct absorption: Simple, if absorption is strong enough
- WMS: More sensitive especially for small signals (near zero baseline)
 - WMS with TDLs improves noise rejection
 - *Normalized WMS, e.g. $2f/1f$ cancels scattering losses!*


3. Absorption Sensor Strategies:

WMS-2f/1f Accounts for Non-Absorption Losses




- Fixed λ WMS-2f/1f
- Ambient H₂O (T=296 K, 60% RH)
- L=29.5 cm, ~6% absorbance)

1392 nm, Partially Blocking Beam



1392 nm, Vibrating Pitch Lens




- Demonstrate normalized WMS-2f/1f in laboratory air
 - 2f/1f unchanged when beam attenuated (e.g., scattering losses)
 - 2f/1f unchanged when optical alignment is spoiled by vibration

WMS-2f/1f signals free of window fouling or particulate scattering

What species/wavelengths can we access?


4. Wavelength Access: A Wide Range of Combustion Species/Applications using Wavelengths in the IR




- Small species such as NO, CO, CO_2 , and H_2O have discrete rotational transitions in the vibrational bands
- Larger molecules, e.g., hydrocarbon fuels, have blended features
- Different strategies used to monitor discrete lines or blended absorption features

4. Wavelength Access:

TDLs Access Visible to Mid-IR



- Allows access to many atoms and molecules
- Visible and telecom TDLs can be fiber-coupled
- TDLs at wavelengths $>$ telecom are emerging rapidly




Now lets consider hardware: lasers and detectors

4. Wavelength Access - Lasers

- Sources – Semiconductor lasers

- Available from the near UV (375 nm) to the far-IR ($\sim 11 \mu\text{m}$)
 - Power: ~ 1 to 500 mW
 - Low power restricts their application to absorption experiments
- Near-IR lasers are compact, rugged, and fiber-coupled
- DFB lasers can be rapidly tuned over several wavenumbers by changing the injection current or laser temperature
 - External cavity diode lasers can be tuned more than 100 cm^{-1}



Diode lasers, near- to extended-near-IR
(\$1000 - \$6000)
Fiber-coupled up to $2.3 \mu\text{m}$



QC lasers, mid-IR ($\sim \$40,000$)



4. Wavelength Access - Detectors


- Detectors – Photodiode/Photovoltaic detectors
 - A photodiode is a semiconductor that generates voltage or current when light is incident on it
 - Like photoconductors, they have a minimum photon energy associated with the bandgap energy of the semiconductor
 - **Source of noise:** Johnson noise (not shot-noise limited)
 - A variation is avalanche photodiode, signal (volts) = constant x intensity

Detector material	λ [μm]
Si	0.2 – 1.1
Ge	0.4 – 1.8
InAs	1.0 – 3.8
InSb	1.0 – 7.0
InSb (77K)	1.0 – 5.6
HgCdTe (77K)	1.0 – 25.0


4. Wavelength Access - Detectors

- Detectors – Select a detector
 - **Criteria:** wavelength, time response, noise, simplicity, cost ...

Wavelength



Bandwidth



- Frequency bandwidth is important for time-resolved measurements
- Bandwidth depends on the detector area, material, temperature, and pre-amplifier gain



4. Wavelength Access - Detectors

- Detectors – Select a detector

- Detector noise is characterized by the detectivity, D^*

$$D^* = \frac{\sqrt{A_{\text{Detector}} \Delta f}}{\text{NEP}}$$

- Δf = bandwidth
 - NEP = noise equivalent power: the amount of the optical power required to equal the magnitude of the detector noise
 - D^* is improved at lower temperatures (cooling)
- The signal-to-noise ratio (SNR) for a measurement dominated by the detector noise can be calculated using:

$$\text{SNR} = \frac{P_{\text{incident}}}{\text{NEP}} = \frac{P_{\text{incident}} D^*}{\sqrt{A_{\text{Detector}} \Delta f}}$$



- Cost and complexity are also important considerations
- Spatially uniform responsivity is also important
 - Smaller and cooled detectors are more uniform



5. Example TDL Applications - Aerospace

1. Subsonic velocity in a laboratory wind tunnel @SU – ambient H_2O
2. Supersonic velocity in a test facility @NASA – H_2O from vitiation
3. Supersonic combustion @UVa – non-uniform T and species on LOS
 - Solutions for non-uniformities (column density and species weighted T)
 - Measurements in H_2 fueled combustor of H_2O column density and $T_{\text{H}_2\text{O}}$
 - Velocity in H_2 fueled combustor
4. Supersonic combustion @ATK – M10 flow capture by model scramjet
5. Scramjet unstart monitor @AFRL

5.1 Subsonic Velocity @ SU: TDLAS Sensors in Wind Tunnel




- TDL sensor of mass flux based on H₂O absorption
- Velocity from Doppler shift of absorption wavelength
- Validate sensor in subsonic wind tunnel w/ ambient H₂O @ Stanford
- 0.5 m/s precision for V in uniform subsonic flow

Next: Test in supersonic-flow facilities at NASA Langley


5.2 Supersonic Velocity@NASA Langley via TDLAS: Direct-Connect Supersonic Combustion Test Facility

Vitiated
inlet air →





DCSCTF: Simulates atmospheric supersonic and hypersonic flight conditions

- $M=2.65$ nozzle with $T_{\text{static}} \sim 990\text{K}$ and $P_{\text{static}} \sim 0.7\text{ atm}$; simulates $M=5$ flight
 - Add optical access to isolator
 - Measure V , T , mass flux



5.2 Supersonic Velocity @NASA via TDLAS

Supersonic test facility at NASA Langley (2009)




Upstream window Translating sensor Downstream window

- Sensor translates to probe vertical and horizontal planes


5.2 Supersonic Velocity @NASA via TDLAS

Supersonic test facility at NASA Langley (2009)



Upstream window Translating sensor Downstream window


Vertical Scan




- Sensor translates to probe vertical and horizontal planes

5.2 Supersonic Velocity @NASA via TDLAS

Supersonic test facility at NASA Langley (2009)



Upstream window Translating sensor Downstream window




- Sensor translates to probe vertical and horizontal planes
- Fast sensor captures start-up transients in V and T

Next: A supersonic combusting flow @ UVa


5.3 Supersonic Combustion @ Uva

CFD Predicts Non-Uniform H₂O Products of Combustion

UVa Combustor




Two-color TDL Sensor
for T and H₂O



NCSU CFD along
TDLAS LOS @ $x/H=6$

$\Phi=0.17$, $y=0.4''$ from injector wall




Path-Average
Value

- Combustion of jet in supersonic cross flow produces non-uniform distribution
- Two-color sensor is scanned transverse to flow to give $f(y)$ data
- But what about variations along LOS(z)?
 - Predicted H₂O density not uniform


5.3 Supersonic Combustion @ Uva

CFD Predicts Non-Uniform H_2O Products of Combustion

UVa Combustor




Two-color TDL Sensor
for T and H_2O



NCSU CFD along
TDLAS LOS @ $x/H=6$

$\Phi=0.17$, $y=0.4''$ from injector wall



Path-Average Values



- But what about variations along $\text{LOS}(z)$?
 - Predicted H_2O density and temperature are both not uniform

What can be done for this problem?


5.3 Supersonic Combustion @ Uva

Non-Uniform T and H₂O Affect Absorption

NCSU CFD along TDLAS
LOS @ $x/H=6$



Simulated absorbance




- Non-uniformity impacts values of absorbance

Solution for non-uniform H₂O: Introduce column density $\sigma_{\text{H}_2\text{O}}$
-Allows direct comparison with CFD


5.3 Supersonic Combustion @ Uva

Non-Uniform T and H₂O Affect Absorption

NCSU CFD along TDLAS
LOS @ $x/H=6$



NCSU CFD along TDLAS
LOS @ $x/H=6$



■ Non-uniformity impacts values of absorbance

Solution for non-uniform H₂O: Introduce column density σ_{H_2O}
-Allows direct comparison with CFD

Solution for non-uniform T: 1. Optimize the line selection
2. Measure species-specific T_{H_2O}
-Extend to multi-species; e.g., T_{CO_2}
-Another metric for new tests of CFD

5.3 Supersonic Combustion @ Uva

Column Density: New Paradigm to Test CFD

- **Problem:** k_v is not constant along LOS, $\tau_v = \frac{I_t}{I_o} = e^{-\int_0^L k_v dx}$
 - **Why:** T , χ_{H_2O} , and ρ are not uniform
 - **New Approach:** Write absorbance α in terms of ρ_{H_2O}

$$\alpha_v = \int_0^L k_v dl = \int_0^L S(T) \Phi_v(T, \rho, \chi) \rho_{H_2O} R_{H_2O} T dl$$

- **Strategy:**



1. Intelligent line selection ($S(T) \propto 1/T$)

Thus, α_v insensitive to non-uniform T


and integrated absorbance (A) becomes

$$A = \int_{-\infty}^{\infty} \alpha_v d\nu = c \int_0^L \rho_{H_2O} dl = c \times \sigma_{H_2O}$$

2. Integrated absorbance yields column density $\sigma_{H_2O} \equiv \int_0^L \rho_{H_2O} dl$
3. Measured σ_{H_2O} new direct test of CFD
4. T_{H_2O} determined by ratio of absorbance on two H_2O lines
also new test of CFD



5.3 Supersonic Combustion at UVa (H₂ Fueled)




- Two-color NIR sensor developed for H₂O and tested at Stanford
 - Confirmed tunnel static temperature and inflow condition for CFD
 - $\sigma_{\text{H}_2\text{O}}$, $T_{\text{H}_2\text{O}}$ & $V_{\text{H}_2\text{O}}$ data to compare with CFD at three measurement planes downstream of fuel injector

Now lets look at an example of the data compared to CFD

5.3 Supersonic Combustion @UVa


TDLAS Results vs CFD

Column Density



Injector side

Temperature




Injector side

- Measured $\sigma_{\text{H}_2\text{O}}$ in best agreement with CFD at $x/H=6$
 - TDLAS suggest greater fuel penetration than CFD predicts
- Measured $T_{\text{H}_2\text{O}}$ in similar agreement with CFD at all x/H
 - Good agreement near injector wall, higher measured T further from wall

Next: Velocity in this combusting flow

5.3 Supersonic Combustion @UVa

TDLAS Velocity via Doppler Shift




5.3 Supersonic Combustion @UVa

TDLAS V_{H_2O} vs CFD for Combusting Flow

Axial Velocity

Core flow velocity
~950m/s

$\phi=0.17$





- H₂O only present as combustion product (thus no V_{H_2O} far from injector wall)
- Large decrease in core-flow velocity near wall, as predicted by CFD
 - Significantly reduced axial velocity due to addition/mixing of fuel
 - Excellent CFD-TDLAS agreement in shape and magnitude

Next: Measurements in an impulse facility

5.4 Supersonic Combustion @ATK

TDLAS in a Model Scramjet @M10



Reflected Shock Tunnel @ ATK GASL
Mach 5-25



- High-speed flow produced by expansion of detonation driven shock wave
- Impulse facility produces a short test time of steady flow (~3ms at M10)
- Tests conducted at M=10 flight condition for H_2 combustion at $\phi \sim 1$

5.4 Supersonic Combustion @ATK

NIR H₂O Sensor for Model Scramjet



Setup:

- 2-color, 3-LOS, fiber-coupled sensor mounted on model

Goals:



- Characterize model capture @M10 in ATK tunnel (TARE run)
- Provide data for H₂ combustion in model to compare with CFD

Results:

- Determined time-resolved T_{H_2O} and σ_{H_2O} to compare with CFD

5.4 Supersonic Combustion @ATK

M10 Combustion: $\phi = 1.03$, $AoA = 7.5^\circ$




Key Results

- Significant combustion product H_2O observed prior to steady test time
 - Evidence of prompt ignition (before steady flow developed)

Next: Can we learn about flow from signal fluctuations?

5.5 Scramjet Unstart Monitor @AFRL

Example: Fluctuations in T Uniformity via TDLAS




AFRL, WPAFB, 2007

- Simultaneous measurements on 4 H_2O lines
 - Two lines for T_{low} and two lines for T_{high}


5.5 Scramjet Unstart Monitor @AFRL:

Sensor Monitors Time-Resolved T_{low} vs T_{high}



Running FFT
of T_{low} & T_{high}

Ratio low-
frequencies



- $T_{\text{low}} \neq T_{\text{high}}$ indicates temperature is not uniform
- Low-frequency fluctuations anticipate inlet unstart
- Fluctuation sensing : A new paradigm for control!



6. TDLAS for Aeropropulsion – Future Trends

- Portable TDL-based sensors useful for T, V, species and mass flux over wide range of conditions, facilities
- Current and future topics:
 - Characterization/maintenance/control of ground-test facilities
 - Emerging applications in flight systems
 - Extension to UV and mid-IR to access new species
 - CO, CO₂, HC's, radicals, NO, NO₂
 - Advanced propulsion concepts: scramjets, detonation engines, rockets



Next Lecture

TDLAS Applications to Energy Conversion

1. Fuel in IC engines – fuel and T
2. H_2O and T in slagging coal gasifier
3. H_2O in NCCC coal gasifier
4. NO and CO in coal-fired boiler exhaust